

STANDARD DIMENSIONS

	PART NUMBER	Material	NOMINAL VOLUME (mL)	OVERFLOW (BRIM) CAPACITY (mL)	A IN (MM)	B IN (MM)	Ø C IN (MM)
FIGURE 1	170-01-0100	PFA	100	165	4.27" (108.5)	2.12" (53.8)	2.19" (55.6)
	170-01-0100-A	PFA	100	165	4.27" (108.5)	2.12" (53.8)	2.19" (55.6)
	170-01-0250	PFA	250	380	5.75" (146.1)	2.12" (53.8)	2.75" (69.9)
	170-01-0250-A	PFA	250	380	5.75" (146.1)	2.12" (53.8)	2.75" (69.9)
	170-01-0500	PFA	500	545	6.33" (160.9)	2.12" (53.8)	3.35" (85.1)
	170-01-0500-A	PFA	500	545	6.33" (160.9)	2.12" (53.8)	3.35" (85.1)
	170-01-1000	PFA	1000	1200	8.60" (218.4)	2.12" (53.8)	3.98" (101.1)
	170-01-1000-A	PFA	1000	1200	8.60" (218.4)	2.12" (53.8)	3.98" (101.1)
	170-01-2000	PFA	2000	2168	9.77" (248.1)	2.12" (53.8)	5.01" (127.3)
	170-01-2000-A	PFA	2000	2240	9.77" (248.1)	2.12" (53.8)	5.01" (127.3)
FIGURE 2	170-01-0050	PFA	50	58	3.31" (84.2)	1.52" (38.7)	1.44" (36.6)
	170-01-0050-A	PFA	50	58	3.31" (84.2)	1.52" (38.7)	1.44" (36.6)




FIGURE 1




FIGURE 2



TITLE: **PURILLEX® PFA BOTTLE DIMENSIONS**

REV. **I**